



Service cart SimplyFlex

Maximum flexibility in your laboratory

The Service cart SimplyFlex offers maximum flexibility and adaptability in the laboratory: users can customise the workstation to suit their needs by replacing or adding individual modules.

HUMAN MOMENTUM. SINCE 1908.

WALDNER



Your advantages

- + **Modular design**
With SimplyFlex, users can customise the workstation themselves by replacing or adding individual modules, ensuring maximum flexibility. Regardless of whether it is an evaluation station, standing workstation, storage option or with an enclosure.
- + **Simple customisation of the layout**
The laboratory layout can be adapted to the respective processes by the user without the need for complex conversion measures, which increases efficiency.
- + **Secure base frame**
The Service cart SimplyFlex offers a base frame with levelling lifting castors for absolute stability.
- + **Shelf close to the floor**
for additional storage space and an organised work surface. The shelf can be supplemented with a counterweight. The service cart can then also be used without a table frame and table top and can accommodate large devices.
- + **Shelves with functional rail**
For attaching laboratory accessories and magnetic LED task lights to provide optimum illumination of the working environment.
- + **Glass splash protection**
An integrated splash guard for additional safety and hygiene in the workplace.
- + **Customised worktop**
Choose from different worktops to customise the worktop to your specific requirements in terms of chemical resistance and cleanability.
- + **Ergonomic adjustment**
Customised ergonomic adjustments for the laboratory staff are possible depending on the configuration level in order to promote health care and a pleasant working environment. Various table frames are available as modules for this purpose, which can be exchanged as required.

Flexible & adaptable

The Service cart SimplyFlex offers maximum flexibility and adaptability in the laboratory: users can customise the workstation by replacing or adding individual modules. The layout can be flexibly adapted to the processes without time-consuming conversion work. With 100% pre-installation, there is no need for assembly on site. Workstations including media supply can be easily rearranged without disturbing neighbouring work areas. Ideal for companies with changing users and processes.

More information:
waldner.de/en/simplyflex

